

MANUAL TILT EXAMINATION CHAIR

Model 1500-CH



RIGHTmed

1500-CH MANUAL TILT EXAMINATION CHAIR

MANUAL TILT CHAIR -

Pneumatically assisted tilting mechanism for easy recline to any position. (Maximum of 40 degrees.)

TILT RELEASE / CONTROL -

One-handed release/control mechanism for easy chair positioning.

POWER LIFT SYSTEM -

Powered by a Swiss made linear actuator for smooth and virtually silent lifting with a 500lb lifting capacity. The chair seat raises to 30.5" (77.50cm) and lowers to 22.5" (57.15cm).

ROTATION / LOCK -

Rotatable to a full 360 degrees, enabling maximum flexibility in patient positioning with locking ability in any position.

FOOTREST -

A large, broad and stable footrest for safe and easy ingress/egress.

ARMRESTS -

Upholstered for comfort, plus individually retractable for patient flexibility.

HEIGHT CONTROL -

Height controls available from slit lamp arm, stand console (when paired with RIGHTmed 1500-ST Stand) or foot switch.

FOOT CONTROL -

A robust foot switch provides raise/lower control for height adjustment.

HEADREST -

Height adjustable, upholstered headrest for maximum patient comfort.

Specifications: base 22.5" (57.15cm) wide x 23" (58.5cm) depth, seat width 18.5" (47cm), footrest 14" (35.5cm) wide x 8.5" (21.6cm), shipping weight: 270lbs (122.5kg) gross 240lbs (109kg) net. With a policy of continual improvement, we reserve the right to make changes to specifications without advance notice.



Distributed by:

Manufactured in Canada by
blue.S4, Inc.
www.blueS4.com

RIGHTmed

5325 Cleveland Street, Suite 303
Virginia Beach, VA 23462

888-224-6012 | www.rightmedical.com