

COMBO UNIT CHAIR AND INSTRUMENT STAND



Model 1000

RIGHTmed

MODEL 1000 COMBO UNIT - CHAIR AND INSTRUMENT STAND

FEATURE, FUNCTION AND VALUE -

The 1000 Combo Unit was specifically designed to meet the needs of customers who require a compact, yet fully functional instrument delivery system for their office. While delivering all of the basic functionality of a standard chair and stand, this unit occupies a small footprint of 36" (91.5cm) wide with no compromise in quality or ease of use.

UPPER (REFRACTOR) ARM -

An ergonomically designed lever mechanism, conforming to natural hand movements, engages a single lock mechanism in three dimensions (forward/backward/vertical). The arm is capable of a vertical range of 12.5" (31.75cm) and is counterbalanced to accommodate instruments up to 20lbs (9.1kg).

THIRD ARM (OPTIONAL) -

With the same ergonomically designed lever mechanism as the refractor arm, the optional third arm has a vertical range of 12.5" (31.75cm) and is counterbalanced to accommodate instruments up to 30lbs (13.6 kg).

COUNTERBALANCED SLIT LAMP ARM -

The lock/release mechanism is conveniently located under the slit lamp arm and allows height adjustment with single handed control. A 10" (26cm) range of travel from 26" (66cm) to 36" (91.5cm) from the floor and is counterbalanced to accommodate instruments up to 60lbs (27.2kg). The switch controlling the chair height is conveniently located on the front of the arm for easy, effortless access.

HALOGEN OVERHEAD LIGHT -

A modern overhead light is a standard feature on the Model 1000 and employs a long-lasting halogen lamp, providing bright, direct light for general lighting or minor procedures. The metal lamp arm construction allows one handed movement with long lasting durability. The lamp has a low voltage (12V) for safety.

Specifications: base width 36" (91.5cm), height to top of pole: 68" (172.75cm), shipping weight: 425lbs (193kg) gross 370lbs (168kg) net. With a policy of continual improvement, we reserve the right to make changes to specifications without advance notice.



Distributed by:

Manufactured in Canada by
blue.S4, Inc.
www.blueS4.com

RIGHTmed

5325 Cleveland Street, Suite 303
Virginia Beach, VA 23462

888-224-6012 | www.rightmedical.com