



INSTRUMENT STAND

Model 1500-ST



RIGHTmed

MODEL 1500-ST INSTRUMENT STAND

The RIGHTmed 1500-ST Instrument Stand is designed for functional and robust performance in the delivery of instruments to patients in the ophthalmic setting. The modern rounded design and all metal construction deliver durable performance with a softer, aesthetically pleasing appearance.

INSTRUMENT CONSOLE -

The console is equipped with three recharging wells and also features a raise/lower switch for chair control, rheostat for overhead light intensity control and main on/off switch.

UPPER (REFRACTOR) ARM -

An ergonomically designed lever mechanism conforms to natural hand movements and engages a single lock mechanism in three dimensions (forward/backward/vertical). The arm is capable of a vertical range of 12.5" (31.75cm) and is counterbalanced to accommodate instruments up to 20lbs (9.1kg).

THIRD ARM (OPTIONAL) -

With the same ergonomically designed lever mechanism as the refractor arm, the optional third arm has a vertical range of 12.5" (31.75cm) and is counterbalanced to accommodate instruments up to 30lbs (13.6 kg).

COUNTERBALANCED SLIT LAMP ARM -

The lock/release mechanism is conveniently located under the slit lamp arm and allows height adjustment with single handed control. A 10" (26cm) range of travel from 26" (66cm) to 36" (91.5cm) from the floor and is counterbalanced to accommodate instruments up to 60lbs (27.2kg). The switch controlling the chair's height is conveniently located on the front of the arm for easy, effortless access (reversible for left handed setup).

HALOGEN OVERHEAD LIGHT -

A modern overhead light is a standard feature on the Model 1500-ST and employs a long-lasting halogen lamp, providing bright, direct light for general lighting or minor procedures. A metal lamp arm construction allows one handed movement with long lasting durability. The lamp is low voltage (12V) for safety.

Specifications: base width 22.5" (57.15cm) wide x 23" (58.5cm) depth, height to top of console 34" (86.4cm), height to top of pole 68" (172.7cm) shipping weight: 270lbs (122.5 Kg) gross 240lbs (109 Kg) net. With a policy of continual improvement, we reserve the right to make changes to specifications without advance notice.



Distributed by:

Manufactured in Canada by
blue.S4, Inc.
www.blueS4.com

RIGHTmed

5325 Cleveland Street, Suite 303
Virginia Beach, VA 23462

888-224-6012 | www.rightmedical.com